

A woman with long, wavy brown hair and a gentle smile is looking towards the camera. She is wearing a white, vertically-ribbed shirt. In her right hand, she holds a clear glass filled with water. In her left hand, she holds a small, light blue, oval-shaped pill. The background is a soft, out-of-focus indoor setting with warm lighting.

***Personalized Medicine to
Improve and Save Lives***

***Improve effectiveness and minimize risk
with MD GeneticPro®***

 **GeneticPro**®
Advanced Genetic Testing • www.MDGeneticPro.com

Impact Patient Health, Safety and Quality-of-Life.

What is Pharmacogenomics?

Pharmacogenomics (PGx) is the relationship between a patient's inherited genetic makeup and their response to pharmaceutical drugs.

Goals of Pharmacogenomics:

Select Responsive Patients

Avoid Adverse Drug Reactions

Maximize Drug Efficacy

When to Pursue Genetic Testing:

- **Prior** to initiation of therapy when using a drug that's **hard to dose**
- **After** the initiation of therapy when a patient's **not achieving** an expected therapeutic goal



Rapid Metabolizer

I can have unexpected responses to certain medications and increased risk for medication interactions.

**No/Sub-optimal response.*



Normal Metabolizer (Extensive Metabolizer)

I have normal enzyme activity and am likely to have expected response to medications.

Consequences of Mismanagement:

Adverse Drug Reactions Over 2 Million/Year













Death 4th in US

Illness 5th in US

Cost of Mortality from ADRs \$177 Billion

Sources: Trescot, et. al. American Society of International Pain Physicians;

FDA Preventable Adverse Drug Reactions: A Focus on Drug Interactions; <http://www.fda.gov/Drugs/DevelopmentApprovalProcess/DevelopmentResources/DrugInteractionsLabeling/ucm110632.htm> Accessed October 27, 2015.

METABOLISM PROFILE	MEDICATION METABOLISM	PARENT DRUG	METABOLITES	POTENTIAL CLINICAL IMPACT
 RAPID	<ul style="list-style-type: none">Significantly increased enzyme activityMetabolizes certain medications at a significantly higher rate than normal	 LOW	 HIGH	<ul style="list-style-type: none">Unexpected therapeutic responseIncreased risk for adverse effectsIncreased risk for medication interactionsNeed for dose adjustmentNeed for medication change
 NORMAL	<ul style="list-style-type: none">Normal enzyme activityMetabolizes certain medications at a normal rate	 NORM	 NORM	<ul style="list-style-type: none">Likely to have expected response to medications
 INTERMEDIATE	<ul style="list-style-type: none">Reduced enzyme activityMetabolizes certain medications at a somewhat lower rate than expected	 INCREASED	 DECREASED	<ul style="list-style-type: none">Minimal effect on therapeutic response and risk for adverse effects
 POOR	<ul style="list-style-type: none">No enzyme activityMetabolizes certain medications at a significantly lower rate than normal	 HIGH	 LOW	<ul style="list-style-type: none">Unexpected therapeutic responseIncreased risk for adverse effectsIncreased risk for medication interactionsNeed for dose adjustmentNeed for medication change

Intermediate Metabolizer

There's minimal effect on my therapeutic response to drugs and risk for adverse effects.

Poor Metabolizer

I have no enzyme activity and can have unexpected therapeutic responses to certain medications, increased risk for adverse effects, and increased risk for medication interactions.

**Toxic to the patient.*

Targeted Therapy



Efficacy and Safety can vary by medication metabolism profile.

MD GeneticPro® Pharmacogenetic Testing :

- Assesses a number of key drug metabolism enzymes¹
- Identifies the underlying genetic characteristics that impact the metabolism of medications commonly prescribed to patients with pain¹
- Simplifies specimen collection with patient-friendly, saliva-based testing⁵

Pain Management Drugs:

- **Only 42% of patients using prescription pain relief medication receive relief.**⁶
- Genetic variations in CYP450 2D6, 3A4, and 3A5 enzymes have been shown to significantly affect responses to various opioids and **alter clinical outcomes, including likelihood of overdose.**⁷⁻⁸
- Prescription opioid pain relievers account for **more overdose deaths than heroin and cocaine combined - 15,000 vs. 4,000.**⁹
- **Approximately 10% of the population are Poor Metabolizers, 7% are Ultra-Rapid Metabolizers and 35% are carriers of a non-functional 2D6 allele, elevating the risk of ADRs when these individuals take multiple drugs.**

Neuroactive Drugs:


- Approximately **50% of all patients do not adequately respond to first line medication** and require lengthy trial and error attempts to identify the appropriate drug.
- Therapeutic benefit typically takes weeks to achieve after stable dosing.
- Testing for genes that are involved in NAD metabolism should be considered when:
 - Starting a patient on a NAS metabolized by CYP2D6, CYP2C19, CYP3A4 or CYP3A5.
 - Pharmacogenetic information exists that suggests concomitant medications may inhibit CYP2D6, CYP2C19, CYP3A4 or CYP3A5 if the patient is taking one or more NAD.
 - Testing for CYP metabolic status can assist physicians in selecting the appropriate NAD choice and dose at the individual patient level.

Improve Outcomes



Clonidogrel (Plavix®):

- **Clonidogrel is the second-highest selling drug in the world** with >90 million patients worldwide.¹⁰
- Up to 30% of patients treated with standard doses of Clonidogrel respond poorly, thus **increasing their risk of recurrent ischemic events**.¹⁰
- 2C19 carriers have more than **50% increased risk of death from cardiovascular causes, MI, or stroke** compared with non-carriers.¹¹
- Patients with variants in 2C19 have a **3.58 times greater risk for major adverse cardiovascular events such as death, heart attack and stroke**.¹¹
- CYP450 2C19 genetic variations are present in **high percentages of ethnic populations**¹²:



Warfarin:

- **Warfarin is the most commonly prescribed anticoagulant in the world.**
- The consequences of incorrect dosage are severe and, in some cases, life threatening.
- **The average bleeding cost associated with Warfarin is \$15,998.**
- Patient response to Warfarin has been shown to be influenced by genetic factors.
- 2C9 & VKORC1 genetic variants explain >50% of the patient variability in Warfarin dose response.
- **Genetic testing could allow Warfarin users to avoid 85,000 serious bleeding events and 17,000 strokes each year...**This could result in a savings of \$1.1 billion annually.
- Studies have shown that Warfarin therapy guided by genetic testing led to **stable anticoagulation in patients 18 days quicker** than those who followed the average dose protocol.

Sample Report



Premier Medical Laboratory Service (PMLS)

6000A Pelham Road, Greenville, SC 29615

Phone: 877-335-2455

Web: www.premedinc.com

Laboratory Director: Peter Zvejnieks, MD

Comprehensive Pharmacogenetic Report Created for: John Doe

Patient: John Doe

DOB: 1/1/1962

Accession #: 2345

Gender: Male

Collection Date: 2/24/2014

Received Date: 2/25/2014



Ordered By: Dr. Smith

Report Generated: 2/27/2014










Patient Medications

Current Medication List: Plavix, Lopressor

Medications Affected by Patient Genetic Results

-  **Plavix (Clopidogrel)** Evidence Level: 2
Significantly Reduced Response to Clopidogrel
Consider alternative therapy. Example of alternative drugs: Prasugrel (contraindicated in TIA/Stroke patients); Ticagrelor; Aspirin; Aspirin plus Dipyridamole.
-  **Lopressor (Metoprolol)** Evidence Level: 1
Normal Sensitivity to Metoprolol
Metoprolol can be prescribed at standard label recommended-dosage and administration. Selection of proper dosage requires individual titration.

Test Details

	Assay	Results	Phenotype	Clinical Consequences
	CYP2C19	*7/*7	Poor Metabolizer	Consistent with a significant deficiency in CYP2C19 activity. Increased risk for side effects or loss of efficacy with drug substrates.
	CYP2C9	*2/*2	Poor Metabolizer	Consistent with a significant deficiency in CYP2C9 activity. Increased risk for side effects or loss of efficacy with drug substrates.
	CYP2D6	*1/*1	Normal Metabolizer	Consistent with a typical CYP2D6 activity. This test did not identify risks for side effects or loss of efficacy with drug substrates.
	CYP3A4	*3/*3	Intermediate Metabolizer	Consistent with a moderate deficiency in CYP3A4 activity. Caution is advised when prescribing narrow therapeutic index drugs. Alternative drugs or dose adjustment may be required if CYP3A inhibitors or inducers are co-prescribed.
	CYP3A5	*3C/*3C	Poor Metabolizer	Consistent with a poor CYP3A5 activity. Caution is advised when prescribing narrow therapeutic index drugs. Alternative drugs or dose adjustment may be required if CYP3A inhibitors or inducers are co-prescribed.
	VKORC1	-1639G>A G/A	Intermediate Warfarin Sensitivity	VKORC1 is the site of action of warfarin. The patient may require a decrease in warfarin dosage.
	Factor II Factor V Leiden	20210G>A GA 1691G>A GG	Increased Thrombosis Risk	The patient's genotypes for Factor V Leiden and Factor II predicts an increased risk for thrombosis. Consider avoiding estrogen-containing preparations. A short course of prophylactic anticoagulation may be considered in high-risk settings such as surgery.
	MTHFR MTHFR	1298A>C AA 677C>T CC	No Increased Risk of Hyperhomocysteinemia	The patient has a normal MTHFR function and no elevation of plasma homocysteine levels is expected. The risk for venous thromboembolism is not increased.
	Apolipoprotein E	ε3/ε4	Increased Risk of Hyperlipidemia/Atherosclerotic Vascular Disease	The patient's has one copy of an abnormal APOE allele, that is associated with an increased risk for hyperlipidemia/atherosclerotic vascular disease.

Sample Report



Premier Medical Laboratory Service (PMLS)

6000A Pelham Road, Greenville, SC 29615

Phone: 877-335-2455

Web: www.premedinc.com

Laboratory Director: Peter Zvejnieks, MD

Potentially Impacted Medications

Cardiovascular Medications		
Standard Precautions	Use With Caution	Consider Alternatives
Carvedilol (Coreg)	Fluvastatin (Lescol)	Clopidogrel (Plavix)
Flecainide (Tambocor)	Warfarin (Coumadin)	
Irbesartan (Avapro)		
Metoprolol (Lopressor)		
Mexiletine (Mexitol)		
Nebivolol (Bystolic)		
Prasugrel (Effient)		
Propafenone (Rythmol)		
Propranolol (Inderal)		
Ticagrelor (Brilinta)		
Timolol (Timoptic)		
Pain Medications		
Standard Precautions	Use With Caution	Consider Alternatives
Codeine (Codeine)	Carisoprodol (Soma)	
Hydrocodone (Vicodin)	Celecoxib (Celebrex)	
Oxycodone (Percocet)	Flurbiprofen (Ansaïd)	
Tramadol (Ultram)	Piroxicam (Feldene)	
Psychotropic Medications		
Standard Precautions	Use With Caution	Consider Alternatives
Aripiprazole (Abilify)	Amitriptyline (Elavil)	
Atomoxetine (Strattera)	Citalopram (Celexa)	
Clozapine (Clozaril)	Clobazam (Onfi)	
Desipramine (Norpramin)	Clomipramine (Anafranil)	
Desvenlafaxine (Pristiq)	Diazepam (Valium)	
Donepezil (Aricept)	Doxepin (Silenor)	
Duloxetine (Cymbalta)	Escitalopram (Lexapro)	
Galantamine (Razadyne)	Imipramine (Tofranil)	
Haloperidol (Haldol)	Phenytoin (Dilantin)	
Iloperidone (Fanapt)	Sertraline (Zoloft)	
Mirtazapine (Remeron)	Tetrabenazine (Xenazine)	
Nortriptyline (Pamelor)	Trimipramine (Surmontil)	
Olanzapine (Zyprexa)		
Paliperidone (Invega)		
Paroxetine (Paxil)		
Perphenazine (Trilafon)		
Pimozide (Orap)		
Risperidone (Risperdal)		
Thioridazine (Mellaril)		
Venlafaxine (Effexor)		
Other Medications		
Standard Precautions	Use With Caution	Consider Alternatives
Darifenacin (Enablex)	Glimepiride (Amaryl)	
Dexlansoprazole (Dexilant)	Glipizide (Glucotrol)	
Esomeprazole (Nexium)	Glyburide (Micronase)	
Fesoterodine (Toviaz)	Tolbutamide (Orinase)	
Lansoprazole (Prevacid)	Voriconazole (Vfend)	
Omeprazole (Prilosec)		
Ondansetron (Zofran)		
Pantoprazole (Protonix)		
Rabeprazole (Aciphex)		
Tacrolimus (Prograf)		
Tamsulosin (Flomax)		
Tolterodine (Detrol)		

Sample Report



Premier Medical Laboratory Service (PMLS)

6000A Pelham Road, Greenville, SC 29615

Phone: 877-335-2455

Web: www.premedinc.com

Laboratory Director: Peter Zvejnieks, MD

Dosing Guidance

- Amitriptyline (Elavil)** Evidence Level: 1
Increased Sensitivity to Amitriptyline (CYP2C19 *7/*7 Poor Metabolizer)
Consider a 50% reduction of recommended amitriptyline starting dose and monitor the plasma concentrations of amitriptyline and nortriptyline to adjust the dose.
- Carisoprodol (Soma)** Evidence Level: 1
Altered Sensitivity to Carisoprodol (CYP2C19 *7/*7 Poor Metabolizer)
Carisoprodol should be used with caution in patients with reduced CYP2C19 activity. Because there is insufficient data to allow calculation of dose adjustment when carisoprodol is prescribed, consider reducing the dose or consider using an alternative therapy.
- Celecoxib (Celebrex)** Evidence Level: 1
High Sensitivity to Celecoxib (CYP2C9 *2/*2 Poor Metabolizer)
Consider starting at half the lowest recommended dose and evaluate response the first week. Be alert to gastrointestinal adverse events. Consider alternative medication for the management of Juvenile Rheumatoid Arthritis.
- Citalopram (Celexa)** Evidence Level: 1
Increased Sensitivity to Citalopram (CYP2C19 *7/*7 Poor Metabolizer)
Consider using citalopram at lower doses and monitor the patient for side effects. Dose escalations over 20 mg/day for CYP2C19 poor metabolizers are not recommended.
- Clobazam (Onfi)** Evidence Level: 1
Increased Sensitivity to Clobazam (CYP2C19 *7/*7 Poor Metabolizer)
In CYP2C19 Poor metabolizers, plasma levels of the active metabolite N-desmethyloclobazam were 5-fold higher than those found in CYP2C19 normal metabolizers. Therefore, the starting dose should be 5 mg/day and dose titration should proceed slowly according to weight. Patients should be titrated initially to 10 mg /day (≤ 30 kg body weight) or 20 mg/day (> 30 kg body weight). If necessary and based upon clinical response, an additional titration to the maximum doses 20 mg/day (≤ 30 kg body weight) or 40 mg/day (> 30 kg body weight) may be started on day 21.
- Clomipramine (Anafranil)** Evidence Level: 1
Increased Sensitivity to Clomipramine (CYP2C19 *7/*7 Poor Metabolizer)
Consider a 50% reduction of recommended clomipramine starting dose and monitor the plasma concentrations of clomipramine and desmethyl-clomipramine to adjust the dose.
- Clopidogrel (Plavix)** Evidence Level: 2
Significantly Reduced Response to Clopidogrel (CYP2C19 *7/*7 Poor Metabolizer)



Premier Medical Laboratory Service (PMLS)

6000A Pelham Road, Greenville, SC 29615


Phone: 877-335-2455

Web: www.premedinc.com

Laboratory Director: Peter Zvejnieks, MD

Patient Reference

PATIENT CARD: CUT OUT AND CARRY IN YOUR WALLET

 For more information: 877-335-2455 www.premedinc.com			APOE	$\epsilon 3/\epsilon 4$	Increased Risk of Hyperlipidemia/Atherosclerotic Vascular Disease
Pharmacogenetic Test Results			Factor II	20210G>A GA	Increased Risk of Thrombosis
Patient: John Doe DOB 1/1/1962			Factor V Leiden	1691G>A GG	
CYP2C19	*7/*7	Poor Metabolizer	MTHFR	677C>T CC	No Increased Risk of Hyperhomocysteinemia
CYP2C9	*2/*2	Poor Metabolizer		1298A>C AA	
CYP2D6	*1/*1	Normal Metabolizer	VKORC1	-1639G>A G/A	Intermediate Warfarin Sensitivity
CYP3A4	*3/*3	Intermediate Metabolizer			
CYP3A5	*3C/*3C	Poor Metabolizer			

Advantages of MD GeneticPro®

Personalized Medicine to Improve & Save Lives

Results

Help predict patients' response to medications.

Assist in clarifying patients' lower than expected clinical response.

Explains patients' higher than expected incidence of adverse effects.

Insights

Manage treatment strategy to maximize safety and efficacy.

Clarify or validate a patient's UDT results.

Individualized Therapy

- More refined diagnostic testing to enhance patient outcomes.
- Identify treatment and dosing options most likely to result in improved quality of life.
- Avoid adverse drug reactions.

Convenience

Patient-friendly saliva or blood-based testing.

Informative and easy-to-read reports.

Assists in informed decision making.

MD GeneticPro® improves clinical outcomes and increases patient safety.



HIPAA Compliant Web Portal Results

Convenient Online Reporting

Our HIPAA Compliant web portal is private, secure and easy for you and your staff to use. Whether you keep physical records, have a paperless office, or have an integrated Electronic Medical Record (EMR) system, our reporting system can seamlessly integrate with your practice.

Premier Medical
LABORATORY SERVICES

Client Outreach Portal

Welcome to Premier Medical Laboratory Web Result Interface. Please enter your credentials to review your result reports.

Please login using the form below

Email Address:
email@example.com

Password:
enter your password here

Login

Secure
Login

Private and HIPAA Compliant

Reports + Order Entry + Search Result Search Help Log

Welcome to Avalon Client Outreach

Last Login: 06/15/2014 09:12:00 AM

Result Search

Patient Fields
First Name: Last Name: Date of Birth (MM/DD/YYYY): MEDICID/Chart #:

Doctor Fields
First Name: Last Name:

Result Fields
Date of Service: From: To: Date Reported: From: To: Specimen Date: From: To: Imported/Created: From: To: Accession: ☐ Abnormal Only ☐ Unprinted Reports ☐ Show Last Login ☐ Inactivated Only

Include orders not yet received by lab: ☐

Convenient
Search

Search by Patient Name or Date of Service, and choose to see only Abnormal results, Unprinted results, or all reports since your last login.

Check	View	Accession #	Doctor (Name) #	Client (Name) #	Patient Name #	Order Date #	Specimen Date #	Status #
<input checked="" type="checkbox"/>	View	10001	SAMPLE PHYSICIAN (C)	SAMPLE HELLERER CENTER (S)	JANE DOE	06/10/2014 10:17:00 AM	06/10/2014 10:17:00 AM	Patient Entered
<input checked="" type="checkbox"/>	View	10002	SAMPLE PHYSICIAN (C)	SAMPLE HELLERER CENTER (S)	JANE DOE	06/10/2014 10:17:00 AM	06/10/2014 10:17:00 AM	Patient Entered
<input checked="" type="checkbox"/>	View	10003	SAMPLE PHYSICIAN (C)	SAMPLE HELLERER CENTER (S)	JANE DOE	06/10/2014 10:17:00 AM	06/10/2014 10:17:00 AM	Patient Entered
<input checked="" type="checkbox"/>	View	10004	SAMPLE PHYSICIAN (C)	SAMPLE HELLERER CENTER (S)	JANE DOE	06/10/2014 10:17:00 AM	06/10/2014 10:17:00 AM	Patient Entered
<input checked="" type="checkbox"/>	View	10005	SAMPLE PHYSICIAN (C)	SAMPLE HELLERER CENTER (S)	JANE DOE	06/10/2014 10:17:00 AM	06/10/2014 10:17:00 AM	Patient Entered
<input checked="" type="checkbox"/>	View	10006	SAMPLE PHYSICIAN (C)	SAMPLE HELLERER CENTER (S)	JANE DOE	06/10/2014 10:17:00 AM	06/10/2014 10:17:00 AM	Patient Entered
<input checked="" type="checkbox"/>	View	10007	SAMPLE PHYSICIAN (C)	SAMPLE HELLERER CENTER (S)	JANE DOE	06/10/2014 10:17:00 AM	06/10/2014 10:17:00 AM	Patient Entered
<input checked="" type="checkbox"/>	View	10008	SAMPLE PHYSICIAN (C)	SAMPLE HELLERER CENTER (S)	JANE DOE	06/10/2014 10:17:00 AM	06/10/2014 10:17:00 AM	Patient Entered
<input checked="" type="checkbox"/>	View	10009	SAMPLE PHYSICIAN (C)	SAMPLE HELLERER CENTER (S)	JANE DOE	06/10/2014 10:17:00 AM	06/10/2014 10:17:00 AM	Patient Entered
<input checked="" type="checkbox"/>	View	10010	SAMPLE PHYSICIAN (C)	SAMPLE HELLERER CENTER (S)	JANE DOE	06/10/2014 10:17:00 AM	06/10/2014 10:17:00 AM	Patient Entered
<input checked="" type="checkbox"/>	View	10011	SAMPLE PHYSICIAN (C)	SAMPLE HELLERER CENTER (S)	JANE DOE	06/10/2014 10:17:00 AM	06/10/2014 10:17:00 AM	Patient Entered
<input checked="" type="checkbox"/>	View	10012	SAMPLE PHYSICIAN (C)	SAMPLE HELLERER CENTER (S)	JANE DOE	06/10/2014 10:17:00 AM	06/10/2014 10:17:00 AM	Patient Entered
<input checked="" type="checkbox"/>	View	10013	SAMPLE PHYSICIAN (C)	SAMPLE HELLERER CENTER (S)	JANE DOE	06/10/2014 10:17:00 AM	06/10/2014 10:17:00 AM	Patient Entered
<input checked="" type="checkbox"/>	View	10014	SAMPLE PHYSICIAN (C)	SAMPLE HELLERER CENTER (S)	JANE DOE	06/10/2014 10:17:00 AM	06/10/2014 10:17:00 AM	Patient Entered
<input checked="" type="checkbox"/>	View	10015	SAMPLE PHYSICIAN (C)	SAMPLE HELLERER CENTER (S)	JANE DOE	06/10/2014 10:17:00 AM	06/10/2014 10:17:00 AM	Patient Entered
<input checked="" type="checkbox"/>	View	10016	SAMPLE PHYSICIAN (C)	SAMPLE HELLERER CENTER (S)	JANE DOE	06/10/2014 10:17:00 AM	06/10/2014 10:17:00 AM	Patient Entered
<input checked="" type="checkbox"/>	View	10017	SAMPLE PHYSICIAN (C)	SAMPLE HELLERER CENTER (S)	JANE DOE	06/10/2014 10:17:00 AM	06/10/2014 10:17:00 AM	Patient Entered
<input checked="" type="checkbox"/>	View	10018	SAMPLE PHYSICIAN (C)	SAMPLE HELLERER CENTER (S)	JANE DOE	06/10/2014 10:17:00 AM	06/10/2014 10:17:00 AM	Patient Entered
<input checked="" type="checkbox"/>	View	10019	SAMPLE PHYSICIAN (C)	SAMPLE HELLERER CENTER (S)	JANE DOE	06/10/2014 10:17:00 AM	06/10/2014 10:17:00 AM	Patient Entered
<input checked="" type="checkbox"/>	View	10020	SAMPLE PHYSICIAN (C)	SAMPLE HELLERER CENTER (S)	JANE DOE	06/10/2014 10:17:00 AM	06/10/2014 10:17:00 AM	Patient Entered

At-A-Glance
Results

Results are easy to view at a glance. Abnormal reports are highlighted in red, normal in gray. A flag alerts you to results already printed. Choose to view individual reports, or several at once.

GeneticPro
Advanced Genetic Testing - www.geneticpro.com

Premier Medical Laboratory Service (PMLS)
6000A Pathway Road, Greenville, SC 29615
Phone: 877-335-2455
Web: www.geneticpro.com

Comprehensive Pharmacogenetic Report Created for: John Doe
Laboratory Director: Peter Zepf, MD

Patient: John Doe DOB: 11/1982
Accession #: 2345 Gender: Male
Collection Date: 2/24/2014 Received Date: 2/25/2014
Ordered By: Dr. Smith Report Generated: 2/27/2014

Patient Medications

Current Medication List: Plavix, Liposol

Medications Affected by Patient Genetic Results

Plavix (Clopidogrel)
Significantly Reduced Response to Clopidogrel
Consider alternative therapy. Examples of alternative drugs: Prasugrel (contraindicated in TIA/stroke patients); Ticagrelor; Aspirin; Aspirin plus Dipyridamol.
Evidence Level: 2

Liposol (Metoprolol)
Normal Sensitivity to Metoprolol
Metoprolol can be prescribed at standard label recommended dosage and administration. Selection of proper dosage requires individualization.
Evidence Level: 1

Test Details

Assay	Results	Phenotype	Clinical Consequences
CYP2C9	*1/*1	Poor Metabolizer	Consistent with a significant deficiency in CYP2C9 activity. Increased risk for side effects or loss of efficacy with drug substrates.
CYP2C9	*2/*2	Poor Metabolizer	Consistent with a significant deficiency in CYP2C9 activity. Increased risk for side effects or loss of efficacy with drug substrates.
CYP2D6	*1/*1	Normal Metabolizer	Consistent with a typical CYP2D6 activity. This test did not identify risks for side effects or loss of efficacy with drug substrates.
CYP3A4	*3/*3	Intermediate Metabolizer	Consistent with a moderate deficiency in CYP3A4 activity. Caution is advised when prescribing narrow therapeutic index drugs. Alternative drug or dose adjustment may be required if CYP3A4 inhibitors or inducers are co-prescribed.
CYP3A5	*3/*3C	Poor Metabolizer	Consistent with a poor CYP3A5 activity. Caution is advised when prescribing narrow therapeutic index drugs. Alternative drug or dose adjustment may be required if CYP3A5 inhibitors or inducers are co-prescribed.
VKORC1	-1639G>A GA	Intermediate Warfarin Sensitivity	VKORC1 is the site of action of warfarin. The patient may require a decrease in warfarin dosage.
Factor II Factor V Leiden	20210GA GA 16910AA GG	Increased Thrombotic Risk	The patient's genotypes for Factor V Leiden and Factor II predicts an increased risk for thrombosis. Consider avoiding estrogen-containing preparations. A short course of prophylaxis (anticoagulation) may be considered in high-risk settings such as surgery.
MTFR MTFR	128AAC AA 877G>T CC	No Increased Risk of Hyperhomocysteinemia	The patient has a normal MTFR function and no elevation of plasma homocysteine levels is expected. The risk for venous thromboembolism is not increased.
Apolipoprotein E	ε3/ε4	Increased Risk of Hyperlipidemia/Atherosclerosis/Vascular Disease	The patient has one copy of an abnormal ApoE allele, that is associated with an increased risk for hyperlipidemia/atherosclerosis/vascular disease.

Genetic Test Results For John Doe Page 1 of 7

Integrated
Reports

Easily save or print PDF reports, or automatically download to your EMR system.



Premier Medical
LABORATORY SERVICES

John Doe

jdoe@premedinc.com

Direct 877-335-2455 ext. XXX

Toll Free 877-335-2455


Fax 877-889-9157

P.O. Box 2009 Fountain Inn, SC 29644 • www.premedinc.com



Personalized Medicine
Save and Improve Lives
Fully Reimbursed by Medicare

1. Argoff CE. Clin J Pain. 2010;26(1 Suppl):S16-S20. 2. Mayo Reference Services Communique. September 2005. <http://www.mayomedicallaboratories.com/media/articles/communique/mc2831-0905.pdf> Accessed June 1, 2012. 3. Crews KR, et al. Clin Pharm Ther. 2012;91(2):321-326. 4. American Medical Association. <http://www.ama-assn.org/resources/doc/genetics/pgx-brochure-2011.pdf>. Accessed June 1, 2012. 5. Koka S, et al. J Am Dent Assoc 2008;139(6):735-740. 6. Peter D Hart Research Assoc, Americans Talk About Pain: A Survey Among Adults Nationwide, Aug. 2003. 7. Janetto, et al. Utilization of Pharmacogenomics and Therapeutic Drug Monitoring for Opioid Management, Oct. 2009. 8. Jin, et al. PGxas Molecular Autopsy for Forensic Toxicology: Genotyping CYP450 3A4*1B and 3A5*3 for 25 Fentanyl cases, Oct. 2005. 9. American Academy of Pain Medicine. 10. Simon, T. et. al. Genetic Determinants of Response to Clopidogrel and Cardiovascular Events. NEJM, Dec, 2009. 11. MedWatch: U.S. Department of Health & Human Services, U.S. F.D.A., <http://www.fda.gov/Safety/MedWatch/default.htm>. 12. Mega, J. et. al, Cytochrome P-450 Polymorphisms and Response to Clopidogrel, NEJM, Dec 2008.



Premier Medical™
LABORATORY SERVICES

www.PreMedInc.com • 877.335.2455 • Fax: 877.889.9157
6000A Pelham Road • Greenville, SC 29615