



Corporate Headquarters: UVP, Inc.
 2066 W. 11th Street, Upland, CA 91786
 Tel: (800)452-6788 * (909)946-3197
 E-Mail: uvp@uvp.com

European Operations: Ultra-Violet Products Ltd.
 Unit 1, Trinity Hall Farm Estate, Nuffield Rd
 Cambridge CB4 1TG UK Tel: +44(0)1223-420022
 E-Mail: uvp@uvp.co.uk

DETERMINING EPROM ERASURE OR STERILIZATION TIME

Sterilization dosage = (UV intensity) x (Exposure time). Therefore,

$$\text{Sterilization time in seconds} = \frac{\text{Sterilization dosage: (W.sec/cm}^2\text{) x 1,000,000}}{\text{Intensity: (\mu W/cm}^2\text{)}}$$

or

$$\frac{\text{Sterilization dosage: (\mu W.sec/cm}^2\text{)}}{\text{Intensity: (\mu W/cm}^2\text{)}}$$

Notice the term “μ” in the numerator of the second equation - this is not in the first equation. By multiplying the numerator in the first equation by 1,000,000, you are converting *watt* into *microwatt*. One microwatt = one watt x 10⁻⁶ (10 to the (-6) power).

Assume you are sterilizing ostrich or emu eggs and the sterilization dosage = 8800 μW.sec/cm²; the intensity of the light source is 1500μW/cm². Note that we are referring to *microwatt* seconds. Since the dosage is already in microwatts, there is no need to multiply the numerator by 1,000,000.

The resulting equation would be:

$$\frac{\text{Sterilization dosage: } 8800\mu\text{W.sec/cm}^2}{\text{UV Intensity: } 1500\mu\text{W/cm}^2} = 5.87 \text{ seconds}$$

Now assume you are erasing EPROMs and the recommended erasure dosage is 15W.sec/cm². Assume the intensity of the light source is 18,000W/cm². The equation would be:

$$\frac{15\text{Wsec/cm}^2 \times 1,000,000}{18,000\mu\text{W/cm}^2 \times 60} \quad \text{You get} \quad \frac{15,000,000\mu\text{Wsec/cm}^2}{18,000\mu\text{W/cm}^2 \times 60} =$$

$$\frac{833.33 \text{ sec.}}{60} = 13.89 \text{ minutes}$$